ESPAÑOL

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ENGLISH

This wheelchair is designed for frequent or occasional use by a single occupant weighing up to 100 kg and can be used indoors and outdoors.

Folding and unfolding the wheelchair: To unfold the wheelchair, open the frame by pulling the armrests, then press down on the side tubes until you are sure that the wheelchair is fully unfolded. Before folding the wheelchair, completely remove the footrests or make sure that they are fixed and in an upright position. To fold the wheelchair, only lift the centre edge of the seat upholstery.

Locking the wheels: Push the brakes forward until they engage in the locked position. Do not attempt to sit down or stand up without first applying the brakes. Always activate both brakes simultaneously.

Adjusting the length of the footrests: The position of the footrests (standard - flat - or raised) can be adjusted to suit the length of the user's legs. The footrest should be adjusted to support the weight of the thighs. The lower edge of the footrest should be at least 7 cm above the floor to allow sufficient clearance when walking on ramps or slopes. To do this, (Maestranza, Giralda, Cathedral), loosen the adjusting screw with a spanner and fold the footrests inwards or outwards to achieve the desired position. Then tighten the screw carefully.

Using the rear support pedals and brakes: Use the rear support pedals to raise the front wheels, for example, when climbing a kerb. To do this, push the rear support pedals down with one foot.

Climbing a kerb: Approach the kerb facing the kerb. The caregiver should use the rear support pedals to raise the front wheels and lower them onto the kerb once the kerb has been raised. Finally, the caregiver should push the wheelchair forward, lifting it slightly to raise the kerb if necessary.

Getting off a kerb: Place the front wheels in front of the kerb edge. The caregiver should use the support pedals to raise the front wheels and tilt the occupant back slightly. While keeping the wheels up, lower the wheelchair slightly towards the kerb.

Warning: When lowering a kerb, the front wheels should be raised to prevent the occupant from falling.

Maintenance:

The user or caregiver should check the safety of the wheelchair at regular intervals.

1. Check the handgrips and non-slip rubbers and ensure that they are solid and secure.
2. For a correct fit, check the wheel locks and make sure that the large wheels are locked when activated.
3. Check all bolts, nuts and joints to ensure that they remain solid.
4. If the wheelchair is fitted with tyres, check that the tyre pressure is correct.
5. Check that the adjusted length of the footrest is correct.
6. Check the rear wheel axle and front fork before using the wheelchair.
7. Check and ensure that all parts are tight and securely locked before using

the wheelchair.

Notice to the user and/or patient

Any serious incident that has occurred in connection with the product should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Indications

- It is designed for the elderly, the patient, the frail and the disabled, as a mobile wheelchair.
- Weight capacity 100 kg.

Tips and precautions

- Do not drive on a bumpy road. Avoid hitting the wheelchair against obstacles.
- When driving the wheelchair on inclined terrain, the user must wear the seat belt. When driving downhill with an angle of inclination of more than 5 degrees, the accompanying person should turn the wheelchair face upwards and drive the wheelchair downhill slowly and carefully. Do not drive the wheelchair downhill with the user upside down or it will cause the wheelchair to tip over and injure the user.
- Do not drop the wheelchair by its own force of inertia onto a step, this could damage the frame or wheels.
- Do not climb up or down steps and do not climb stairs without the help of an assistant. Do not lift the wheelchair with the user seated in it, otherwise you will damage the wheelchair and cause injury to the user or depth.
- Do not lean your body out of the wheelchair to prevent the wheelchair from tipping over.
- Do not drive the wheelchair on escalators.
- Do not step on the footrest when getting in and out of the wheelchair. Do not stand on the footrest. Otherwise it could tip over.
- Do not tilt the wheelchair or turn the wheelchair suddenly. Do not brake the wheelchair suddenly.
- The parking brake is only for parking the wheelchair. Do not use it to slow down the wheelchair or to brake the wheelchair in motion.
- Do not force the chassis or push handle to fold the wheelchair. To fold the wheelchair, lift the seat of the wheelchair from the centre.
- Do not lean forward in the wheelchair unless both feet are on the floor. This is especially important for users with heavy casts with raised foot rests.

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GUANGDONG DAYANG MEDICAL TECHNOLOGY CO., LTD.
SUNGO Europe B.V.
Imported by/Importado por: Grupo R. Queraltó S.A.